



Gems[®]

**DESIGNED BY
EMBRYOLOGISTS
WHO USE IT¹**

Gems[®]: A complete embryo medium suite, including the third generation of the world's first sequential media^{1,2}





Gems®

For embryos cultured outside of the uterus, *in vivo* conditions must be recreated to stimulate embryo development.³ An **optimal culture environment** is vital to support the overall performance of the ART lab and is essential to minimise stress to gametes and embryos.³⁻⁵

In the 1990s, Genea IVF clinic in Australia (formerly Sydney IVF) developed one of the world's first sequential human embryo culture media suites.¹⁻²

INTRODUCING Gems[®]

YOUR COMPLETE RANGE OF IVF AND VITRIFICATION MEDIA

Gems[®] is the third generation of this culture medium suite ^{1, 2} developed originally by Sydney IVF/Genea IVF clinic and further refined by Genea Biomedx in-house scientists. It has been conceived to support your needs at every stage of the ART process – from gamete collection right through to embryo transfer and cryopreservation.



OPTIMISED FOR EVERY STAGE:

Gems[®] Medium Suite includes gamete and embryo handling buffers and solutions for embryo culture and cryopreservation. ^{6, 7, 8}



ROBUST AND EFFECTIVE MEDIA:

A wide variety of amino acids, nutrients, vitamins, balanced salts and unique antioxidants to support embryo development in vitro. ⁶



LONG HERITAGE:

Over 30-year history and experience in culture media development, aiming for optimum fertility outcomes. ^{1, 2}



ESTABLISHED SAFETY AND EFFICACY:

Gems[®] media has been in clinical use at Genea Fertility clinics in Australia since 2013, helping more than 200,000 embryos to be cultured and more than 16,000 babies born. ⁹



CONSISTENTLY HIGH QUALITY BATCH-TO-BATCH:

Manufactured in state-of-the-art ISO13485-certified facilities in Sydney, Australia, undergoing extensive Quality Control testing before release. ⁶




USER-FRIENDLY:

Gamete and Culture media are ready-to-use with clear labeling with distinctive colours and icons to aid ease of recognition in the lab.

DEVELOPED BY USERS FOR USERS

OVER 30 YEARS DEVELOPING CULTURE MEDIA

Experts from the fields of IVF and chemistry have contributed to the development of Gems® media suite.





2011



Gems® clinical trial
ran at Genea IVF
clinics in 2011-2013

2013



3rd Generation media
- Gems - in full clinical
use at Genea IVF clinics.

2016



Gems® media
manufacturing started
at ISO13485 - certified
facility in Sydney, Australia.

2016



Gems® media suite
obtains CE Mark.

2017



Majority of Gems®
media suite obtains
FDA clearance.*

* Sperm solutions pending FDA clearance.

OPTIMISED FOR EVERY STEP

Gems® adapts to your laboratory practice and workload demands with a complete range of IVF and vitrification media with differing bottle sizes.

GAMETE HANDLING AND PREPARATION



Oocyte Retrieval Buffer

ORB-50 (50mL)
ORB-20 (20mL)



- HEPES-buffered, supplemented with gentamycin (0.01 mg/mL), not containing human serum albumin.



Sperm Wash Gradient Set

SWG-45P-50 (50mL)
SWG-90P-50 (50mL) (not for individual sale)



- Used for isolating motile sperm and removing seminal plasma from semen in preparation for insemination.
- HEPES-buffered, supplemented with gentamicin (0.01 mg/mL), human serum albumin (10 mg/mL) and coated silica.



Sperm Medium

SPM-50 (50mL)
SPM-20 (20mL)



- Used for washing sperm for IUI, IVF, ICSI and diagnostics.
- Bicarbonate-buffered, supplemented with gentamicin (0.01 mg/mL) and human serum albumin (10 mg/mL).



Sperm Buffer

SPM-50 (50mL)
SPM-20 (20mL)



- Used for washing sperm for IUI, IVF, ICSI and diagnostics.
- HEPES-buffered, supplemented with gentamicin (0.01 mg/mL) and human serum albumin (10 mg/mL).

VITRIFICATION & WARMING SOLUTIONS



Vitrification Set

VIT-SOL1-20 (20mL)
VIT-SOL2-20 (20mL)
VIT-SOL3-10 (10mL)
(not for individual sale)

- Used for vitrifying embryos.
- HEPES buffered, supplemented with Ethylene Glycol, Dimethyl Sulfoxide (DMSO), trehalose and human serum albumin (20mg/mL).



Warming Set

VRM-SOL1-20 (20mL)
WRM-SOL2-10 (10mL)
WRM-SOL3-10 (10mL)
(not for individual sale)

- Used for warming vitrified embryos.
- HEPES buffered, supplemented with trehalose and human serum albumin (20mg/mL).

OTHERS



VitBase

VBS-20 (20 mL)

- Used for maintaining embryos for a short period of time in a non-gassed environment.
- HEPES buffered, supplemented with human serum albumin (20mg/mL).

OPTIMISED FOR EVERY STEP

CULTURE MEDIA

Each Gems® culture medium has a tailored composition for optimal support of embryos, including amino acids, vitamins, nutrients, balanced salts and a pool of antioxidants, which are key for normal metabolism of gametes and embryos. ^{2,3} Sequential media formulations are developed to adapt to developing embryos' needs, containing EDTA and low levels of glucose and essential amino acids for early stages, and no EDTA and higher levels of glucose and essential amino acids for later stages. Geri Medium has been developed to support embryos throughout their in vitro development period for uninterrupted culture requiring no media change or replenishment ¹⁰. It contains medium levels of EDTA, glucose and essential amino acids. All Gems® culture media are optimised for use in a low oxygen (5%) environment. ^{2,3}



Fertilisation Medium

FEM-50 (50mL)
FEM-20 (20mL)

- Used for providing a suitable environment for fertilisation.
- Bicarbonate-buffered, supplemented with gentamicin (0.011mg/mL), human serum albumin (5 mg/mL) and high concentration of glucose.



Cleavage Medium

CLM-50 (50mL)
CLM-20 (20mL)

- Used for culturing embryos from fertilisation to Day 3.
- Bicarbonate-buffered, supplemented with gentamicin (0.011mg/mL), human serum albumin (5 mg/mL), antioxidants and low concentration of glucose.



Blastocyst Medium

BLM-50 (50mL)
BLM-20 (20mL)

- Used for culturing embryos from Day 3 to blastocyst stage.
- Bicarbonate-buffered, supplemented with gentamicin (0.011mg/mL), human serum albumin (5 mg/mL), antioxidants and high concentration of glucose.



Geri® Medium


ONE-50 (50mL)
ONE-20 (20mL)



Used for culturing embryos uninterrupted from fertilisation to blastocyst stage.

- Bicarbonate-buffered, supplemented with gentamicin (0.011mg/mL), human serum albumin (5 mg/mL), antioxidants and medium concentration of glucose.
- Contains Glutamine in dipeptide form as Alanyl-Glutamine to prevent harmful ammonium build-up in medium.
- Helps simplify workflow and processes by eliminating the need for medium change or use of additional culture dishes during culture.

Whether your lab employs sequential or uninterrupted continuous embryo culture, with Gems[®] you have the flexibility to choose your preferred culture method.



CONSTITUENTS & DISTINGUISHING COMPONENTS



Gems® media contains key components to support embryo growth.

Gems® medium suite products are buffered with bicarbonate or HEPES. They include gentamicin for microbial protection and varying levels of human serum albumin (HSA) as a protein source, depending on the product and its role in the IVF cycle.

There are many benefits to HSA in culture media¹²⁻¹³ :

- Aiding embryo metabolism.
- Protecting against toxins.
- Assisting in pH buffering.
- Acting as a colloid osmotic regulator.
- Preventing embryos and gametes sticking to the device used to handle and culture them



Gems® culture media are advanced formulations including antioxidants to support embryo development.

- **L-carnitine:** An antioxidant that protects DNA and cell membranes from free radical damage and from apoptosis.¹⁸

L-carnitine supplementation has also been shown to have a positive effect in cryosurvival rates in bovine and porcine embryos.^{16,17}

L-carnitine supplementation has been recently correlated with higher number of good quality embryos and improved clinical outcomes in human embryos.¹⁸


- **Cobalmin (vitaminB₁₂):** Important factor in methionine synthase reaction, critical for maintaining normal DNA methylation¹⁹ processes.
- **Folic acid (vitaminB₉):** An antioxidant that supports cell cleavage and growth. It is required for synthesis, repair and protection of DNA during methylation.^{20,21,22}
- **Ascorbic acid (vitamin C):** An antioxidant that protects DNA and cell membranes from oxidative damage and apoptosis, reduces DNA fragmentation and abnormal gene expression, and when co-incubated with embryos, improves blastocyst development ²³

THE PROOF IS IN THE RESULTS

Gems® media have been in clinical use at Genea IVF clinics since 2013 across all cycle types and patient ages. Implementation of a fully undisturbed culture system using Geri® Medium alongside Geri® incubator resulted in an increase in blastocyst quality, embryo utilisation and pregnancy rates⁹.

Over 200,000 embryos have been cultured and over 16,000 babies been born with the help of Gems media⁹.

MAXIMISE THE POTENTIAL OF UNDISTURBED INCUBATION



* Statistically significant differences ($p < 0.01$).


**Genea modified Gardner grading.

*** Zygotes utilised prior to Day 5 were removed from the total # of embryos cultured because their potential fate was unknown.

A statistically significant increase ($p < 0.001$) in live births/ongoing pregnancies was seen when Genea IVF clinics changed culture system from benchtop incubator and sequential media to Geri® timelapse incubator and Geri® Medium. This shows that not only embryo development is improved by the optimised micro- environment of uninterrupted culture, but the practice translates to better clinical outcomes.²⁴

GENEA IVF CLINICS' RESULTS OF ALL CYCLES, INCLUDING CUMULATIVE OUTCOMES FROM CYCLES FOLLOWED FOR UP TO 2 YEARS"

Live births/ongoing pregnancies per cycle *



- Geri® & Geri® Medium (Aug 2014 - Dec 2015; Total # embryos cultured 22,928)
- MINC® & Sequential Gems® (Aug 2016 - Dec 2017; Total # embryos cultured 19,964)

* Transition rate from foetal heart pregnancy approx 94% regardless of the culture system.

EASE OF USE

Pre-supplemented with human serum albumin (HSA) as a protein source ^{2, 13},
Gems supports standardised laboratory practices.

Ready-to-use solutions
to aid standardisation

Different bottle sizes
adapt for different
size labs and needs.

Clear symbols and
colour coded labels
aid recognition.



EXTENSIVE QUALITY AND REGULATORY CONTROL

Critical media ingredients and materials are screened before used for manufacturing. Finished products undergo applicable QC tests using analytical methods validated according to international pharmacopoeia and adhering to stringent acceptance criteria before product release to ensure their safety.²⁵



All Gems[®] media suite products have CE Mark, and most of them have FDA clearance*. Most media suite products are available in >40 countries.³⁵

* Sperm solutions pending FDA clearance.

REFERENCES

- 1- Mortimer D. Human blastocyst development media. Hum Reprod. 2001 Dec;16(12):2725. doi: 10.1093/humrep/16.12.2725.
- 2- Mortimer & Mortimer. Essential features in media development for spermatozoa, oocytes, and embryos. In 'Culture Media, Solutions, and Systems in Human ART', Ed P Quinn. 2014 (pp. 47–67). Cambridge University Press.
- 3- Swain et al. Optimizing the culture environment and embryo manipulation to help maintain embryo developmental potential. Fertil Steril. 2016 Mar;105(3):571-587. doi: 10.1016/j.fertnstert.2016.01.035.
- 4- Testori S. Embryo culture media update (2015). HFEA Scientific and Clinical Advances Advisory Committee (SCAAC), SCAAC (10/51)01.https://ffqlive.blob.core.windows.net/umbraco-website/1831/embryo_culture_media_update.pdf.
- 5- Will et al. Biological pH buffers in IVF: help or hindrance to success. J Assist Reprod Genet. 2011 Aug;28(8):711-24. doi: 10.1007/s10815-011-9582-0.
- 6- Genea Biomedx Data on File (QFRM190 Gems® Clinical Evaluation Report)
- 7- QFRM40 Gems® instructions for use
- 8- QFRM866 Gems® instructions for use (US)
- 9- Genea Biomedx Data on File (QRTV318 Human Embryo Culture in Geri®)
- 10- Rusden et al. Impact of media replenishment for day six extended blastocyst culture in single-step culture. FSAANZ Conference 14-19.9 2024, Perth, WA, Australia
- 11- Hashimoto et al. Medium without ammonium accumulation supports the developmental competence of human embryos. J Reprod Dev. 2008 Oct;54(5):370-4. doi: 10.1262/jrd.20012.
- 12- Blake et al. Protein supplementation of human IVF culture media. J Assist Reprod Genet. 2002 Mar;19(3):137-43. doi: 10.1023/a:1014788821965.
- 13- Laverge et al. Prospective randomized study comparing human serum albumin with fetal cord serum as protein supplement in culture medium for in-vitro fertilization. Hum Reprod. 1997 Oct;12(10):2263-6. doi: 10.1093/humrep/12.10.2263.
- 14- Abdelrazik et al. L-carnitine decreases DNA damage and improves the in vitro blastocyst development rate in mouse embryos. Fertil Steril. 2009 Feb;91(2):589-96. doi: 10.1016/j.fertnstert.2007.11.067.
- 15- Mansour et al. L-carnitine supplementation reduces oocyte cytoskeleton damage and embryo apoptosis induced by incubation in peritoneal fluid from patients with endometriosis. Fertil Steril. 2009 May;91(5 Suppl):2079-86. doi: 10.1016/j.fertnstert.2008.02.097.
- 16- Lowe et al. Supplementation of culture medium with L-carnitine improves the development and cryotolerance of in vitro-produced porcine embryos. Reprod Fertil Dev. 2017 Nov;29(12):2357-2366. doi: 10.1071/RD16442.
- 17- Takahashi et al. Supplementation of culture medium with L-carnitine improves development and cryotolerance of bovine embryos produced in vitro. Reprod Fertil Dev. 2013;25(4):589-99. doi: 10.1071/RD11262.
- 18- Kim et al. Effects and pregnancy outcomes of L-carnitine supplementation in culture media for human embryo development from in vitro fertilization. J Obstet Gynaecol Res. 2018 Nov;44(11):2059-2066. doi: 10.1111/jog.13763.
- 19- Halczuk et al. Vitamin B12-Multifaceted In Vivo Functions and In Vitro Applications. Nutrients. 2023 Jun 13;15(12):2734. doi: 10.3390/nu15122734.
- 20- Crider et al. Folate and DNA methylation: a review of molecular mechanisms and the evidence for folate's role. Adv Nutr. 2012 Jan;3(1):21-38. doi: 10.3945/an.111.000992.
- 21- Antony. In utero physiology: role of folic acid in nutrient delivery and fetal development. Am J Clin Nutr. 2007 Feb;85(2):598S-603S. doi: 10.1093/ajcn/85.2.598S.
- 22- O'Neill. Endogenous folic acid is essential for normal development of preimplantation embryos. Hum Reprod. 1998 May;13(5):1312-6. doi: 10.1093/humrep/13.5.1312.
- 23- Wang et al. Vitamin C and vitamin E supplementation reduce oxidative stress-induced embryo toxicity and improve the blastocyst development rate. Fertil Steril. 2002 Dec;78(6):1272-7. doi: 10.1016/s0015-0282(02)04236-x.
- 24- Peura et al. Is the biggest impact on clinical IVF outcomes obtained by implementation of continuous media, time-lapse incubator or both? 2020 Fertil Steril, 114:3, e129
- 25- Ph Eur: Pharmacopoeia European; BP: British Pharmacopoeia; USP: US Pharmacopoeia; CP: Chinese Pharmacopoeia; IP: Indian Pharmacopoeia; NF: National Formulary (US)
- 26- Current BP (XVI A); Current Ph Eur (2.6.1); Current USP/NF <71> (External Lab)
- 27- Sterility: External Assay, Current BP (XVI A); Current Ph Eur (2.6.1); Current USP/NF <71>
- 28- pH: Inhouse validated method as per current BP, Ph Eur, USP/NF
- 29- Osmolality: Inhouse validated method as per current BP, Ph Eur, USP/NF
- 30- MEA: Inhouse validated 1-cell MEA method incl. positive and negative controls
- 31- SPTT: Inhouse validated method incl. two independent test samples and negative control
- 32- Gentamicin: Inhouse validated Gentamicin ELISA kit Assay
- 33- Endotoxin: Inhouse validated LAL (limulus ameobocyte lysate) test as per Current BP, Ph Eur, USP/NF
- 34- HSA: Inhouse validated BCA Protein Assay
- 35- Genea Biomedx Data on File (QFRM89 Master Product List)

For healthcare professionals only. Please refer to the instructions for use.

Gems® Culture medium suite complies with applicable medical device standards.

Gems® is manufactured by Genea Biomedx.

Gems Culture medium suite complies with applicable medical device standards.

www.geneabiomedx.com



Product complies with applicable
European Union (EU) regulations

