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1. Swain, Jason E., *et al.* "Optimizing the culture environment and embryo manipulation to help maintain embryo developmental potential". *Fertility and sterility* 105.3 (2016): 571-587.
2. Zhang, J. Q., *et al.* Reduction in exposure of human embryos outside the incubator enhances embryo quality and blastulation rate. *Reproductive Biomedicine Online* 2010; 20(4):510-515.
3. Swain, Jason E. "Decisions for the IVF laboratory: comparative analysis of embryo culture incubators". *Reproductive Biomedicine Online* 28.5 (2014): 535-547.
4. QFRM422 Geri® User Manual.
5. Bonteke, S., *et al.* Low oxygen concentrations for embryo culture in assisted reproductive technologies. *Cochrane Database of Systematic Reviews* 2012, Issue 7.
6. Kirkegaard, K., *et al.* "Effect of oxygen concentration on human embryo development evaluated by time-lapse monitoring". *Fertility and sterility* 99.3 (2013): 738-744.

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1. Swain, J. E., et al. "Optimizing the culture environment and embryo manipulation to help maintain embryo developmental potential". Fertility and sterility 105.3 (2016): 571-587.

2. Zhang, J. Q., et al. "Reduction in exposure of human embryos outside the incubator enhances embryo quality and blastulation rate". Reproductive Biomedicine Online 20(4):510-515.

3. Swain, Jason E. "Decisions for the IVF laboratory: comparative analysis of embryo culture incubators". Reproductive Biomedicine Online 28.5 (2014): 535-547.


4. QFRM422 Geri® User Manual.

5. Takenaka, M., et al. "Effects of light on development of mammalian zygotes". Proceedings of the National Academy of Sciences 104.36 (2007): 14289-14293.

6. Schultz, R. M. "Of light and mouse embryos: less is more". Proceedings of the National Academy of Sciences 104.37 (2007): 14547-14548.

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8. OFRM584 Geri® Clinical Evaluation Report.

9. Swain, J. E., *et al.* "Media osmolality changes over 7 days following culture in a non-humidified benchtop incubator". Fertility and Sterility 106.3 (2016): e362.

10. Swain, J. E., *et al.* "Microdrop preparation factors influence culture-media osmolality, which can impair mouse embryo preimplantation development". Reprod Biomed Online. 2012; 24:142-147.

11. Napier, C, *et al.* "Abstract "Is humidification of incubators necessary for uninterrupted culture?" ACE 9th Biennial Conference 2014.

12. Albert, C, *et al.* "The effect of high humidity culture conditions over embryo development: a continuous embryo monitoring assessment". Reproductive BioMedicine Online 37 (2018): e15-e16.

13. de los Santos JM, *et al.*, "The effect of high humidity culture conditions over embryo development: a continuous embryo monitoring assessment", ESHRE 2019 abstract.

14. Perez Albala, S. *et al.*, "Embryo culture conditions under high humidity significantly enhances blastocysts formation and quality according to an automatic time-lapse algorithm", ESHRE 2019 abstract.

15. Fawzy, M., *et al.* "Humid versus dry incubator: A prospective, randomized, controlled trial". Fertility and sterility 108.2 (2017): 277-283.

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
1. Kovacs, P. (2014) "Embryo selection: the role of time-lapse monitoring". Reprod Biol Endocrinol 12 (1): 124-135.

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

1. Data on file at Genea Biomedx. QRTV318_Human Embryo Culture in Geri®.

2. QFRM422 Geri® User Manual.

3. Swain, J. E., et al. "Optimizing the culture environment and embryo manipulation to help maintain embryo developmental potential". Fertility and Sterility 105.3 (2016): 571-587.

4. Zhang, J. Q., et al. "Reduction in exposure of human embryos outside the incubator enhances embryo quality and blastulation rate". Reproductive Biomedicine Online 20(4):510-515.

5. Swain, J. E., et al. "Decisions for the IVF laboratory: comparative analysis of embryo culture incubators". Reproductive Biomedicine Online 28.5 (2014): 535-547.

6. Bontekoe, S., et al. "Low oxygen concentrations for embryo culture in assisted reproductive technologies". The Cochrane Library (2012).

7. Kovacs, P., et al. "Embryo selection: the role of time-lapse monitoring". Reprod Biol Endocrinol 12 (1) (2014): 124-135.

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
9. De los Santos JM, et al., "The effect of high humidity culture conditions over embryo development: a continuous embryo monitoring assessment", ESHRE 2019 abstract.

10. Fawzy, M., et al. "Humid versus dry incubator: A prospective, randomized, controlled trial". Fertility and sterility 108.2 (2017): 277-283..

11. Perez Albala, S., et al., "Embryo culture conditions under high humidity significantly enhances blastocysts formation and quality according to an automatic time-lapse algorithm", ESHRE 2019 abstract.

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
1. QFRM794 Geri® Connect & Geri® Assess User Manual.
2. Resumen de las pruebas con QRTF285_03 Geri 6.0 y GC&GA 2.0.
3. ALPHA scientists in Reproduction medicine and ESHRE SIG of Embryology (2011). Istanbul consensus workshop on embryo assessment: Proceedings of an expert meeting. Human Reproduction 26 (6): 1270-1283
4. Ciray, H.N., *et al.* (2014). Proposed guidelines on the nomenclature and annotation of dynamic human embryo monitoring by a time-lapse user group. Human Reproduction 29 (12): 2650-2660.
5. Meseguer, M., *et al.* "Using artificial intelligence (AI) and time-lapse to improve human blastocyst morphology evaluation". Human Reproduction (2018): 125-126.
6. Basile, N., *et al.* (2015). The use of morphokinetic1. ics as a predictor of implantation: a multicentric study to define and validate an algorithm for embryo selection. Human Reproduction 30 (2): 276-283.

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

1. 358_1053_01_Instructions for Use Eeva® System 3.1.

2. QFRM794 Geri® Connect and Geri® Assess User Manual.

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4. Adamson, G. D., et al. Improved implantation rates of day 3 embryo transfers with the use of an automated time-lapse-enabled test to aid in embryo selection. *Fertility and sterility* 105.2 (2016): 369-375.

5. Aparicio-Ruiz, B., et al. Automatic time-lapse instrument is superior to single-point morphology observation for selecting viable embryos: retrospective study in oocyte donation. *Fertility and Sterility* 106.6 (2016): 1379-1385.

6. Rocafort, E., et al. "Euploid embryos selected by an automated time-lapse system have superior SET outcomes than selected solely by conventional morphology assessment". *Journal of assisted reproduction and genetics* 35.9 (2018): 1573-1583.

7. Behr, B., et al. Non-invasive technology combining time-lapse Imaging and statistical modeling: Bringing automation into the lab to improve blastocyst selection. *ASRM*. 2015.

8. Pool, TB. et al., "Computer vs. human manual assessment: are we ready to start an automated era in embryology laboratory?", *Fertility and Sterility* 105.2 (2016): e37

9. Aparicio-Ruiz, B., et al. "Clinical validation of an automatic time-lapse algorithm classification system for blastocyst selection", Abstract accepted at ESHRE 2019

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