Genea BIOMEDX



Genie Solution

The Genie Solution provides comprehensive genetic analysis of pre-implantation embryo biopsy samples including the detection of aneuploidy, structural rearrangements and linkage analysis to detect the inheritance of single gene disorders. End-to-end solution where, sample libraires are prepared with Genie Kits, sequenced on the Genie Sequencer followed by data analysis and reporting on the Genie Data Analysis Software.





For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE-A SOLUTION

Genie-A is a solution for detecting CNV, ploidy, Sibling QC and ROH (regions of homozygosity) in preimplantation embryo biopsy samples. We offer two approaches: Genie-A Fast and Genie-A Advanced.

KEY FEATURES

COMPREHENSIVE DETECTION.

A single solution to detect CNV, Ploidy, Sibling QC and ROH (regions of homozygosity).

DEFAULT SNP ANALYSIS.

Genie-A uses genome-wide SNP information to detect Ploidy, Sibling QC and ROH.

HIGH RESOLUTION.

≥4Mb CNVs (unknown), and ≥1 Mb inherited CNVs (known).

HIGH THROUGHPUT AND SCALABILITY.

5M reads/samples, 1x100 bp reads, Up to 48, 96 samples/run.

Genie-A Kits	Analysis Methods	Reporting Content
Genie-A Fast	• CBS Algorithm • Log Likelihood Ratio (LLR) • SNPs BAF Analysis	• Aneuploidy, Ploidy (Triploidy, Haploidy) • ≥4Mb CNVs in unknown regions • ≥1Mb known inherited CNVs • Mosaicism ≥30% C ≥10Mb • Sibling QC • Whole chromosome level ROH
Genie-A Advanced		

For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE-A: TWO APPROACHES



GENIE-A FAST:

Fast workflow, combining WGA and library preparation to reduce cost and time.




GENIE-A ADVANCED:


Separate WGA (MDA) followed by a patented library preparation method to support wider genome coverage.

TWO VERSIONS OF GENIE-A KITS WITH DIFFERENT SAMPLE PREPARATION METHODS FOR GREATER FLEXIBILITY, BUT THE SAME OUTCOME.


ACCURATE DETECTION OF CNV AND PLOIDY



del(7)(p22.3p22.1)(4.5Mb), +10



Triploid (5M reads)



UPD (5M reads)

For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE-PLUS SOLUTION

Genie-Plus Solution is a genome-wide sequencing solution designed to detect the inheritance of monogenic disorders, structural rearrangements, in addition to CNV, ploidy, UPD, Sibling QC and regions of homozygosity (ROH) in preimplantation embryo biopsy samples.

KEY FEATURES

COMPREHENSIVE DETECTION.

A single solution to detect the inheritance of monogenic disorders, structural rearrangements, in addition to CNV, Ploidy, UPD, sibling QC and ROH.

GENOME-WIDE COVERAGE.

WGA by MDA and a patented library preparation method for better whole-genome coverage.

HIGH RESOLUTION.

≥4Mb CNVs (unknown), and ≥1 Mb inherited CNVs (known).

HIGH THROUGHPUT AND SCALABILITY.

40M reads/samples, 2x100 bp reads, Up to 24, 48 samples/run.




For Research Use Only (RUO). Not for use in diagnostic procedures.


COMPLETE GENETIC ANALYSIS

The Genie-Plus solution performs comprehensive genetic analysis of embryo biopsy samples and provides results including, haplotype chart, CNV chart, BAF chart, ROH chart, LLR chart, Z-Score chart and Kinship analysis.

Haplotype chart




CNV chart




BAF chart




ROH chart




Kinship analysis



LLR chart



Z-Score chart



Genie Kits	Analysis Methods	Reporting Content
Genie-Plus	• CBS Algorithm • Log Likelihood Ratio (LLR) • SNPs BAF Analysis	• Aneuploidy, Ploidy (Triploidy, Haploidy) • ≥4Mb CNVs in unknown regions • ≥1Mb known inherited CNVs • Mosaicism ≥30% C ≥10Mb • Sibling QC • Whole chromosome level ROH
	• Linkage Analysis • SNP analysis • CNV Analysis	• > 4000 monogenic Disorders • Balanced translocations, Robertsonian translocations, inversion

For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE SECURE DATA ANALYSIS SYSTEM

Highest security level and privacy during data analysis, data generated by the Genie Sequencer is transferred to the secure local server, through a LAN connection for data analysis on the Genie Data Analysis Software.

- **MAXIMUM SECURITY:**
  All data analysis happens locally on a secure local server, eliminating cloud risks and ensuring full data protection.

- **100% LOCAL, 100% PRIVATE:**
  Your data never leaves your infrastructure, ensuring complete control and compliance with strict security standards.

WITH GENIE, YOU GET CUTTING-EDGE TECHNOLOGY WITH UNMATCHED SECURITY



LOCAL LAB

For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE SEQUENCER

The Genie Sequencer supports the sequencing of sample libraries prepared through the Genie Kits.

- Compact benchtop sequencer designed for ease of use and integration.
- High-throughput sequencer generates 500M reads/run, offering the ability to be used for various applications, including carrier screening, whole exome sequencing, and whole genome sequencing.
- Patterned flow cell facilitates highly accurate sequencing with Q30 >85%, decreased duplicates, and reduced index hopping.




Genie Solution	No. of Samples/run	Reads/sample	Read Length	Application
Genie-A Fast	Up to 48 & 96	5M	1 x 100 bp	CNV, Ploidy, Sibling QC, ROH
Genie-A Advanced	Up to 48 & 96	5M	1 x 100 bp	CNV, Ploidy, Sibling QC, ROH
Genie-Plus	Up to 24 & 48	40M	2 x 100 bp	Linkage analysis, SR, CNV, Ploidy, Sibling QC, ROH



For Research Use Only (RUO). Not for use in diagnostic procedures.


GENIE: SINGLE SOLUTION

A single patented workflow for the detection of CNV, ploidy, Sibling QC, ROH (Genie-A Advanced), structural rearrangements and linkage-based analysis to detect the inheritance of single gene disorders (Genie-Plus).



GENIE: END TO END SOLUTION

Genie provides an end-to-end solution for the genetic analysis of embryo biopsy samples. Sample libraires are prepared using the Genie kits, sequenced on the Genie Sequencer, followed by data analysis and reporting on the Genie Data Analysis Software.



For Research Use Only (RUO). Not for use in diagnostic procedures.

GENIE KIT CATALOGUE

For Research Use Only (RUO). Not for use in diagnostic procedures.

Catalogue Number	Product Name	Packing List	Quantity
GEN-GAF-96-1	**Genie-A Fast** for Genie Sequencer (96 samples/run)	Genie-A Sample Prep Kit (96 reactions)	1
		Genie Sequencing kit	1
		Genie Sequencing Flow Cell	1
GEN-GAF-48-2	**Genie-A Fast** for Genie Sequencer (48 samples/run x2)	Genie-A Sample Prep Kit (96 reactions)	1
		Genie Sequencing kit	2
		Genie Sequencing Flow Cell	2
GEN-GAA-96-1	**Genie-A Advanced** for Genie Sequencer (96 samples/run)	Genie WGA kit (48 reactions)	2
		Genie-Plus Sample Prep Kit (50 reactions)	2
		Genie-Plus Index Primers (96 indexes)	1
		Genie Sequencing kit	1
		Genie Sequencing Flow Cell	1
GEN-GAA-48-1	**Genie-A Advanced** for Genie Sequencer (48 samples/run)	Genie WGA kit (48 reactions)	1
		Genie-Plus Sample Prep Kit (50 reactions)	1
		Genie-Plus Index Primers (96 indexes)	1
		Genie Sequencing kit	1
		Genie Sequencing Flow Cell	1
GEN-GAP-24-2	**Genie-Plus** for Genie Sequencer (24 samples/run x2)	Genie WGA kit (48 reactions)	1
		Genie-Plus Sample Prep Kit (50 reactions)	1
		Genie-Plus Index Primers (96 indexes)	1
		Genie Sequencing kit	2
		Genie Sequencing Flow Cell	2
GEN-GAP-24-4	**Genie-Plus** for Genie Sequencer (24 samples/run x4)	Genie WGA kit (48 reactions)	2
		Genie-Plus Sample Prep Kit (50 reactions)	2
		Genie-Plus Index Primers (96 indexes)	1
		Genie Sequencing kit	4
		Genie Sequencing Flow Cell	4
GEN-WGA-01*	**Genie WGA kit**	Genie WGA kit (48 reactions)	1

**Common to both Genie-A Advanced and Genie-Plus*

For detailed catalogue options to support DNBSEQ-G400, please contact your Genea representative.

For Research Use Only (RUO). Not for use in diagnostic procedures.

Please refer to the instructions for use.

For further information, please contact your Service Representative or visit:

www.geneabiomedx.com

MRK-0014842-02

For Research Use Only (RUO). Not for use in diagnostic procedures.